

ME-380S

Multifunction medical chair for blood donor & dialysis dialysis chair&blood donor chair



Our **blood donation bed** require very little effort to use. We've designed the chair to work with a patient's body, not against it. And your healthcare staff will appreciate our design as well. Less stooping, easier access for cleaning and locking front brakes are among the many features.

- Functions of chair**
- 1, Height is adjustable from 560~800mm
 - 2, Back section is adjustable from 0° ~ 70°
 - 3, Knee section is adjustable from -45° ~ 0°
 - 4, can be a bed

5, Multiple positions available

6, Three German motors or Darnmak Linak hospital motors

7, loading capacity 250kg

8, PU standard for American CA117 or British BS5852 that is resistant to blood, sweat, fire, urine, spit or scraping.

Optional accessories :

1. Instrument tray
2. Infusion stand

1 Humanization design

Modern design ,high quality materials

Modern detailed finishing of the upholstery as well as mechanical and electrical parts underline the high standard of the comfort-Line. Continuous adjustment to all sitting and lying positions, adjustable armrest, utmost comfort for your patients .It makes the chairs suitable for dialysis centers, blood bank, blood donor centers and all kind of luxury hospital rooms.

Muti-position adjustment

Option of choosing three motors for maximum patient comfort and reduction in staff work load. Backrest, seatrest, legrest, height are fully adjustable .

2 Saving space

Only 1.2Sq.Mtr.space required ,saving about 40% space than the tradition medical bed .

3 Dust-proof & easily disinfect design

Good ABS covering for protecting the chair not corrosion ,waterproof ,dustproof ,easily clean & disinfect & maintain .

4 High quality products

Taking the aero seat as the prototype improvement ,upholstery with polyurethane ,luxurious. comfortable and durability . And warm & light color makes patient enjoy the treatment .