



Shanghai Medeco Industry Co.,Ltd.

No.222 Meidan Road Shanghai, China.201908

Tel: 0086-21-32509306 Fax: 0086-21-51686997

Mob: 0086-13761409730 Email:info@medeco-medical.com

www.medeco-medical.com

<http://mermaid-medical.en.made-in-china.com>

Operating lamp Operating Table Hospital Bed
Manufacturer OEM service Medical Solutions



Hydraulic manual Operating Table 1088 new type



Hydraulic manual **Operating Table1088** is widely used in chest

Abdominal surgery, ENT, gynaecology and obstetrics, urology and orthopedics etc.

Hydraulic lift by foot pedal, head operated movements.

The base and column cover all made by premium 304 stainless steel.

Tabletop is made of Composite laminate for X-ray, makes a High definition image.

SPECIFICATION

Table height: 760mm~1010mm

Length and width: 2010mmx500mm

Reversed Trendelenburg/Trendelenburg: 22°

Lateral Tilt (Left/Right): 20°

Head rest(up/down): 60°/70°

Back rest (up/down): 75°/15°



Shanghai Medeco Industry Co.,Ltd.

No.222 Meidan Road Shanghai, China.201908

Tel: 0086-21-32509306 Fax: 0086-21-51686997

Mob: 0086-13761409730 Email:info@medeco-medical.com

www.medeco-medical.com

<http://mermaid-medical.en.made-in-china.com>

Operating lamp Operating Table Hospital Bed
Manufacturer OEM service Medical Solutions



Foot Board (up/down&outward): 15° /85° / 90°

Kidney Bridge up (Optional):120mm

Standard Accessories

1. Pads :1 Set
2. Leg support:1 Pair
3. Arm support:1Pair
4. shoulder support:1 Pair
5. Body support:1 Pair
6. Anesthetic Screen :1 pc
7. Footboard:1Pair

