



Shanghai Medeco Industry Co.,Ltd.

No.222 Meidan Road Shanghai, China.201908

Tel: 0086-21-32509306 Fax: 0086-21-51686997

Mob: 0086-13761409730 Email:info@medeco-medical.com

www.medeco-medical.com

<http://mermaid-medical.en.made-in-china.com>

Operating lamp Operating Table Hospital Bed
Manufacturer OEM service Medical Solutions



Electric Operating Table DT-12F new type



Operating Table DT-12F New type as an advanced automatic operating table, can meet multi-purpose requirements of surgery operations, long horizontal sliding, suitable for C-arm and X-Ray.

It adopted micro touch remote control; any movements can be adjusted by it. Flexible adjustment on head board, back board and seat plate, built-in kidney bridge.

High automation, low noise, high reliability.

Key parts imported from abroad, can be regarded as an ideal electric table.



Shanghai Medeco Industry Co.,Ltd.

No.222 Meidan Road Shanghai, China.201908

Tel: 0086-21-32509306 Fax: 0086-21-51686997

Mob: 0086-13761409730 Email:info@medeco-medical.com

www.medeco-medical.com

<http://mermaid-medical.en.made-in-china.com>

Operating lamp Operating Table Hospital Bed
Manufacturer OEM service Medical Solutions



SPECIFICATION

Table height: 700mm~1000mm

Length and width: 2010mmx500mm

Reversed Trendelenburg/Trendelenburg: 25°

Lateral Tilt (Left/Right): 15°

Head rest(up/down): 45° /90°

Back rest (up/down): 75° /20°

Foot Board (up/down&outward): 15° /90° / 90°

Kidney Bridge up :110mm

Horizontal sliding:300mm

Power:220V/50Hz,1.0Kw

