



**Shanghai Medeco Industry Co.,Ltd.**

No.222 Meidan Road Shanghai, China.201908

Tel: 0086-21-32509306 Fax: 0086-21-51686997

Mob: 0086-13761409730 Email:info@medeco-medical.com

[www.medeco-medical.com](http://www.medeco-medical.com)

<http://mermaid-medical.en.made-in-china.com>

Operating lamp Operating Table Hospital Bed  
Manufacturer OEM service Medical Solutions



## **ELECTRIC GENERAL GYNECOLOGY OPERATING TABLE 3004 (B) ECOC005**



Medical equipment multi-purpose gynecological table is designed for delivery, abortion and gynecological examination. Its characterized by easy operation, safely and reliability, economical use, it is composed of bed surface, bedstead and bed base, the surface includes back section, seat section and leg section, the up-down movement of back section and the forward and back tilt of bed surface are controlled by the operation hand wheel, which makes doctor and patients feel comfortable.

### **SPECIFICATION**

Power supply: 220V/50HZ

Power: 150VA

Table height: 680mm~880mm

Length and width: 1780mmx600mm

Reversed Trendelenburg: 80°

Trendelenburg: 20°



**Shanghai Medeco Industry Co.,Ltd.**

No.222 Meidan Road Shanghai, China.201908

Tel: 0086-21-32509306 Fax: 0086-21-51686997

Mob: 0086-13761409730 Email:info@medeco-medical.com

[www.medeco-medical.com](http://www.medeco-medical.com)

<http://mermaid-medical.en.made-in-china.com>

Operating lamp Operating Table Hospital Bed  
Manufacturer OEM service Medical Solutions



Upward folding angle of back rest:45°

Downward folding angle of back rest: 10°

Rated load weight: 135KG

#### Standard Accessories

1. Pads
2. Leg bracket (Foot placing holder)
3. Hand placing holder
4. Pull Handle
5. Feculence(Filtration screen )