



Shanghai Medeco Industry Co.,Ltd.

No.222 Meidan Road Shanghai, China.201908

Tel: 0086-21-32509306 Fax: 0086-21-51686997

Mob: 0086-13761409730 Email:info@medeco-medical.com

www.medeco-medical.com

<http://mermaid-medical.en.made-in-china.com>

Operating lamp Operating Table Hospital Bed
Manufacturer OEM service Medical Solutions



ELECTRIC GYNECOLOGY EXAMINATION TABLE 3004 (ECOC006)



Power supply: 220V/50HZ

Input Power: 90VA

Environment temperature: 5°C~40°C

Relative humidity: ≤80%

Atmospheric pressure: 860hPa~1060hPa

Back board size: 830mmx550mm

Buttock board size: 400mmx500mm

Table board adjustable range (Except Mattress): Minimum 600mm±20mm

Up and down distance:≥280mm

Upward folding angle of back rest: 70°

Buttock board adjusting range:

Reversed Trendelenburg≥ 8°

Trendelenburg≥ 15°

Standard Accessories



Shanghai Medeco Industry Co.,Ltd.

No.222 Meidan Road Shanghai, China.201908

Tel: 0086-21-32509306 Fax: 0086-21-51686997

Mob: 0086-13761409730 Email:info@medeco-medical.com

www.medeco-medical.com

<http://mermaid-medical.en.made-in-china.com>

Operating lamp Operating Table Hospital Bed
Manufacturer OEM service Medical Solutions



1. Legs bracket 2 sets
 2. Arm support 1set
 3. Reel of paper 1 roll
 4. Feculence basin 1 sets
 5. Remote control 1pc
 6. Foot Switch 1pc
- Option: Dental stool