



Shanghai Medeco Industry Co.,Ltd.

No.222 Meidan Road Shanghai, China.201908

Tel: 0086-21-32509306 Fax: 0086-21-51686997

Mob: 0086-13761409730 Email:info@medeco-medical.com

www.medeco-medical.com

<http://mermaid-medical.en.made-in-china.com>

Operating lamp Operating Table Hospital Bed
Manufacturer OEM service Medical Solutions



Side-control mechanical operating table 3001 Series



Model 3001: medical equipment Side-control mechanical operating table (table surface, base and pump shell with painting).

Model 3001A: medical equipment Side-control mechanical operating table (base and pump shell with stainless steel).

Model 3001B: Medical equipment Side-control mechanical operating table (table surface, base and pump shell with stainless steel).

3001 series comprehensive operating table is necessary for various hospital to operate, which is designed for operation on head, neck, abdomen, chest, the four limbs and others in general surgery. Using the current latest oil pump for hydraulic elevation, and the other operations are controlled manually.

SPECIFICATION

Table height: 800mm~1045mm

Length and width: 2100mmx480mm

Reversed Trendelenburg: 35°



Shanghai Medeco Industry Co.,Ltd.

No.222 Meidan Road Shanghai, China.201908
Tel: 0086-21-32509306 Fax: 0086-21-51686997
Mob: 0086-13761409730 Email:info@medeco-medical.com
www.medeco-medical.com
<http://mermaid-medical.en.made-in-china.com>

Operating lamp Operating Table Hospital Bed
Manufacturer OEM service Medical Solutions



Trendelenburg: 20°

Left inclining angle: 22°

Right inclining angel: 22°

Upward folding angle of back rest: 75°

Downward folding angle of back rest: 22°

Upward folding angle of head rest: 90° (Dismountable, extendable)

Downward folding angle of head rest: 90°

Downward folding angle of leg rest: 90° (Dismountable)

Oil pump range: 240mm

Option: Kidney bridge range=120mm

Standard Accessories

1. Pads
2. Leg bracket
3. Arm rest
4. Shoulder rest
5. Body support
6. Anesthetic Screen
7. Body restraint strap