



**Shanghai Medeco Industry Co.,Ltd.**

No.222 Meidan Road Shanghai, China.201908

Tel: 0086-21-32509306 Fax: 0086-21-51686997

Mob: 0086-13761409730 Email:info@medeco-medical.com

[www.medeco-medical.com](http://www.medeco-medical.com)

<http://mermaid-medical.en.made-in-china.com>

Operating lamp Operating Table Hospital Bed  
Manufacturer OEM service Medical Solutions



## Mechanical hydraulic operating table 3008H



### DESCRIPTION

**Back rest** and Leg rest controlled by air cylinder imported from Germany.

**Soft PVC**, antistatic, latex free, radiolucent pads are also available in a new anti-decubitus material, which meets the safety standard, the surface is anti-bacteria, easily clean, according with disinfection needs.

**Mechanical hydraulic operating table 3008H** pedestal and other main accessories, Rail, Base Shell and Column shell are made of Stainless Steel, Antistatic. Entirely made in stainless steel,

**The surface of Mechanical hydraulic operating table 3008H** is plated by anti-static cover for easy cleaning and preventing from scratching and collision, Nice sculpt supply enough impulsion resistance and stability, easy to be clean and disinfect without dead angle leftover. Special Lower height 650mm be suited for Ophthalmic-operation

### SPECIFICATION

Table height: 650mm~890mm



**Shanghai Medeco Industry Co.,Ltd.**

No.222 Meidan Road Shanghai, China.201908  
Tel: 0086-21-32509306 Fax: 0086-21-51686997  
Mob: 0086-13761409730 Email:info@medeco-medical.com  
[www.medeco-medical.com](http://www.medeco-medical.com)  
<http://mermaid-medical.en.made-in-china.com>

Operating lamp Operating Table Hospital Bed  
Manufacturer OEM service Medical Solutions



Length and width: 2000mm x 480mm  
Reversed Trendelenburg: 28 °  
Trendelenburg: 28 °  
Left inclining angle: 20 °  
Right inclining angel: 20 °  
Upward folding angle of back rest: 75 °  
Downward folding angle of back rest:10 °  
Upward folding angle of head rest: 35 °  
Downward folding angle of head rest: 90 °  
Downward folding angle of leg rest:90 °

### **Standard Accessories**

1. Pads
2. Leg bracket
3. Arm rest
4. shoulder rest
5. Body support
6. Anesthetic Screen
7. Body restraint strap