



Shanghai Medeco Industry Co.,Ltd.

No.222 Meidan Road Shanghai, China.201908

Tel: 0086-21-32509306 Fax: 0086-21-51686997

Mob: 0086-13761409730 Email:info@medeco-medical.com

www.medeco-medical.com

<http://mermaid-medical.en.made-in-china.com>

Operating lamp Operating Table Hospital Bed
Manufacturer OEM service Medical Solutions



A-5 Three function manual bed with ABS headboards

size: 2140mm(L)x940mm(W)x530mm/720mm(H)



Main features:

1. Epoxy coated bed frame, with Collapsible aluminium alloy side rails, and detachable ABS engineering plastic headboards ,
2. three function, which can be operated by manual crank, the details is as follows:
(1)back rest: $0^{\circ} \sim 70^{\circ} \pm 5^{\circ}$ (2)knee rest: $0^{\circ} \sim 35^{\circ} \pm 5^{\circ}$ (3) high-low:530mm~720mm
3. protection strip in bed lateral side & Anti-bumper in four corners, location for IV rod in each corner
4. Four dia-125mm individual castors with brakes