



**Shanghai Medeco Industry Co.,Ltd.**

No.222 Meidan Road Shanghai, China.201908

Tel: 0086-21-32509306 Fax: 0086-21-51686997

Mob: 0086-13761409730 Email:info@medeco-medical.com

[www.medeco-medical.com](http://www.medeco-medical.com)

<http://mermaid-medical.en.made-in-china.com>

Operating lamp Operating Table Hospital Bed  
Manufacturer OEM service Medical Solutions



## **A-10 Three-function Manual Bed**



size: 2140mm(L)x940mm(W)x530mm/720mm(H)

### **Main features:**

1. Epoxy coated bed frame, high quality ABS engineer plastic head/foot board & side rails;
2. three function, which can be operated by manual crank, the details is as follows:  
(1)back rest:0 °70 °±5 ° (2)knee rest:0 °35 °±5 ° (3) high-low:530mm~720mm
3. Anti-bumper in four corners, location for IV rod in each corner
4. Four dia-125mm castors, with central locking system;