



Shanghai Medeco Industry Co.,Ltd.

No.222 Meidan Road Shanghai, China.201908

Tel: 0086-21-32509306 Fax: 0086-21-51686997

Mob: 0086-13761409730 Email:info@medeco-medical.com

www.medeco-medical.com

<http://mermaid-medical.en.made-in-china.com>

Operating lamp Operating Table Hospital Bed
Manufacturer OEM service Medical Solutions



ECOG028-CBII obstetric bed



ECOG028-CBII obstetric bed are designed for various stages of care for patients .They can be used in labor, delivery and postpartum recovery period .The chair posture can be easily converted to provide adequate comfort for the lying-in women .The lying-in women don't need to be moved from one bed to another thanks to the beds' perfect structure.

Features:

Imported LINAK motor Various actions are controlled by remote controller.

Super mute big foot wheel which is easily moved and equipped with brake device.

Equipped with stainless steel sewage basin .

The mattress adopts high density memory sponge and is once formed,with no seam.

Standard Accessories :

A pair of leg support frame,

a stainless steel sewage basin,

a stand for transfusion,

a remote control,

a power wire .