



Shanghai Medeco Industry Co.,Ltd.

No.222 Meidan Road Shanghai, China.201908

Tel: 0086-21-32509306 Fax: 0086-21-51686997

Mob: 0086-13761409730 Email:info@medeco-medical.com

www.medeco-medical.com

<http://mermaid-medical.en.made-in-china.com>

Operating lamp Operating Table Hospital Bed
Manufacturer OEM service Medical Solutions



ECOG030-2E integrated obstetric bed



ECOG030-2E integrated obstetric bed is suitable for conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc.by medical unit.

Features:

Imported LINAK motor Various actions are controlled by remote controller.

The table, base seat and upright column covering are made of high quality 304 stainless steel, which is anti-corrosive and can be easily cleaned.

The mattress adopts high density memory sponge and is once formed,with no seam.

Whole stainless steel hidden auxiliary table with flexible extension and retraction.

Equipped with stainless steel sewage basin and other accessories.

Foot treadle brake device which is easy for moving obstetric table.

Standard Accessories

A pair of arm resting frame,

a pair of leg support frame,



Shanghai Medeco Industry Co.,Ltd.

No.222 Meidan Road Shanghai, China.201908

Tel: 0086-21-32509306 Fax: 0086-21-51686997

Mob: 0086-13761409730 Email:info@medeco-medical.com

www.medeco-medical.com

<http://mermaid-medical.en.made-in-china.com>

Operating lamp Operating Table Hospital Bed
Manufacturer OEM service Medical Solutions



a pair of armrest,
a stainless steel sewage basin,
a hand-holding remote controller.