



**Shanghai Medeco Industry Co.,Ltd.**

No.222 Meidan Road Shanghai, China.201908

Tel: 0086-21-32509306 Fax: 0086-21-51686997

Mob: 0086-13761409730 Email:info@medeco-medical.com

[www.medeco-medical.com](http://www.medeco-medical.com)

<http://mermaid-medical.en.made-in-china.com>

Operating lamp Operating Table Hospital Bed  
Manufacturer OEM service Medical Solutions



## **Electric Hydraulic Operating Table ECOG 006II**

Electric Hydraulic Operating Table ECOG 006II is suitable for various operations conducted by medical unit.



### **Features:**

The bed surface adopts imported carbon plastic plate, upright column is side-lying design fulfills the needs of C type arm whole body examination, photography.

The base seat and upright column covering are made of high quality 304 stainless steel which can be easily cleaned and sterilized.

Both sides of the bed rail for the high quality 304 stainless steel.

The mattress is wrapped with imported artificial leather. High-density memory sponge is formed in one piece, with no seam.

Various actions are controlled according to electric hydraulic working principle and by remote controller.



**Shanghai Medeco Industry Co.,Ltd.**

No.222 Meidan Road Shanghai, China.201908

Tel: 0086-21-32509306 Fax: 0086-21-51686997

Mob: 0086-13761409730 Email:info@medeco-medical.com

[www.medeco-medical.com](http://www.medeco-medical.com)

<http://mermaid-medical.en.made-in-china.com>

Operating lamp Operating Table Hospital Bed  
Manufacturer OEM service Medical Solutions



## **SPECIFICATION:**

Table height: 750mm~950mm

Length and width: 2000mmx520mm

Reversed Trendelenburg/Trendelenburg: 18°/30°

Lateral Tilt (Left/Right): 18°

Head rest(up/down): 50°/90°

Back rest (up/down): 50°/15°

FootBoard downstarts:90°

Kidney bridge range: 120mm

Power: 220V/50~60Hz; 600W