



**Shanghai Medeco Industry Co.,Ltd.**

No.222 Meidan Road Shanghai, China.201908

Tel: 0086-21-32509306 Fax: 0086-21-51686997

Mob: 0086-13761409730 Email:info@medeco-medical.com

[www.medeco-medical.com](http://www.medeco-medical.com)

<http://mermaid-medical.en.made-in-china.com>

Operating lamp Operating Table Hospital Bed  
Manufacturer OEM service Medical Solutions



## MANUAL UNIVERSAL OPERATING TABLE ECOG015-JYC

ECOG015-JYC operation table is suitable for various operations conducted by medical unit.



### Features:

Various actions are controlled by hydraulic pneumatic system and mechanical driven system.

Head plate can be exchanged with leg plate which can fulfill the needs of C type arm complete examination, X-ray photography.

Imported Y type sealing ring is durable.

The table surface adopts imported carbon plastic plate, the mattress adopts high density memory sponge and is once formed, with no seam.

Quality stainless steel 304 is used for base plate, column shield and rails on both sides of the table, easy for cleaning and sterilization.

Foot treadle brake device which is easily for moving operating table. **Standard Accessories:** A pair of shoulder support frame, a pair of waist support frame, a pair of arm resting frame, a pair of leg support frame, a screen frame, a pair of pedal foot board, a waistband.

### SPECIFICATION

Table height: 700mm~1050mm



**Shanghai Medeco Industry Co.,Ltd.**

No.222 Meidan Road Shanghai, China.201908

Tel: 0086-21-32509306 Fax: 0086-21-51686997

Mob: 0086-13761409730 Email:info@medeco-medical.com

[www.medeco-medical.com](http://www.medeco-medical.com)

<http://mermaid-medical.en.made-in-china.com>

Operating lamp Operating Table Hospital Bed  
Manufacturer OEM service Medical Solutions



Length and width: 520mmx2050mm

Reversed Trendelenburg/Trendelenburg: 15 °

Lateral Tilt (Left/Right): 15 °

Head rest(up/down): 30 90 °

Leg board Downwards:90 °