



Shanghai Medeco Industry Co.,Ltd.

No.222 Meidan Road Shanghai, China.201908

Tel: 0086-21-32509306 Fax: 0086-21-51686997

Mob: 0086-13761409730 Email:info@medeco-medical.com

www.medeco-medical.com

<http://mermaid-medical.en.made-in-china.com>

Operating lamp Operating Table Hospital Bed
Manufacturer OEM service Medical Solutions



MANUAL UNIVERSAL OPERATING TABLE ECOG014-JYD



ECOG014-JYD operation tables are suitable for various operations conducted by medical unit.

Features:

Various actions are controlled by hydraulic pneumatic system and mechanical driven system.

The table can be rotated 360° on the condition that the base is fixed.

Imported Y type sealing ring and is durable.

The table surface adopts imported carbon plastic plate which is easily for conducting X-ray examination, photography.

Both sides of the bed rail for the high quality 304 stainless steel.

JY D Table support is formed with alloy aluminum by forging process and will never get rusty.

The mattress is wrapped with imported artificial leather. High-density memory sponge is formed in one piece, with no seam.

JY D The base seat and upright column covering are made of high quality 304 stainless steel which can be easily cleaned and sterilized.

Foot treadle brake device which is easily for moving operating table. Standard Accessories: A pair of



Shanghai Medeco Industry Co.,Ltd.

No.222 Meidan Road Shanghai, China.201908

Tel: 0086-21-32509306 Fax: 0086-21-51686997

Mob: 0086-13761409730 Email:info@medeco-medical.com

www.medeco-medical.com

<http://mermaid-medical.en.made-in-china.com>

Operating lamp Operating Table Hospital Bed
Manufacturer OEM service Medical Solutions



shoulder support frame, a pair of waist support frame, a pair of arm resting frame, a pair of leg support frame, a screen frame, a pair of pedal foot board, a waistband.

SPECIFICATION

Table height:700mm~1050mm

Length and width: 2050mmx520mm

Horizontal rotary:360 °

Reversed Trendelenburg/Trendelenburg: 15 °/20 °

Backboard upwards/Downwards:60 °/10 °

Lateral Tilt (Left/Right): 25 °

Head rest(up/down): 30 °/90 °

Leg board Downwards:90 °

Leg of the board split:180 °