



**Shanghai Medeco Industry Co.,Ltd.**

No.222 Meidan Road Shanghai, China.201908

Tel: 0086-21-32509306 Fax: 0086-21-51686997

Mob: 0086-13761409730 Email:info@medeco-medical.com

[www.medeco-medical.com](http://www.medeco-medical.com)

<http://mermaid-medical.en.made-in-china.com>

Operating lamp Operating Table Hospital Bed  
Manufacturer OEM service Medical Solutions



## ME-360 Full-digital Portable Ultrasound Scanner



### Main Features

1. High-precision full-digital imaging technology, progressive scanning without flickering, delicate and clear images
2. **12.1 inch LED** high definition monitor. Only **6.5kg** Easy for taking, portable, lighter
3. Professional embedded ultrasonic platform
4. high density probe, broadband and frequency conversion technology, Increase the image quality greatly
5. Backlight silica gel keyboard, Photoelectricity track control
6. Intellectualized TGC gain control of the whole paragraph, precisely adjustment image density
7. Wide applications: Abdominal, OB/GYN, Urology, Cardiology, Small parts examination etc.
8. **USB port**: storable and readable

### Technical Specification



**Shanghai Medeco Industry Co.,Ltd.**

No.222 Meidan Road Shanghai, China.201908

Tel: 0086-21-32509306 Fax: 0086-21-51686997

Mob: 0086-13761409730 Email:info@medeco-medical.com

[www.medeco-medical.com](http://www.medeco-medical.com)

<http://mermaid-medical.en.made-in-china.com>

Operating lamp Operating Table Hospital Bed  
Manufacturer OEM service Medical Solutions



Scanning mode: convex / linear / micro-convex

Probes: **96 elements** convex probe,trans-vaginal probe, trans-rectal probe, high frequency linear probe

Operation panel: user-friendly, convenient and flexible with back-lit silica gel keyboard & track ball operation

Monitor: 12.1 inch high resolution LED

Display mode: B, B+B, B+M, M, 4B

Detecting depth:  $\geq 242$ mm

USB port: storable and readable

Image storage

Body marks

Probe interface: 2

Magnification: 16kinds

Cine loop: 256 frames, realized in automatic successive playback and manual single step forward, backward playback

Standard configuration

12.1 inch high resolution LED

3.5MHz Multi-freq Convex Probe(2.5MHz-5.0MHz)

USB port

Options

6.5MHz Trans-vaginal probe(5.5MHz-7.5MHz)

7.5MHz High frequency linear probe(6.5MHz-8.5MHz)

7.5 MHz Multi-freq rectal Linear Probe

Trolley, Printer

- Previous: DW-500 12.1 inch LED Full-digital Laptop Ultrasound Scanner
- Next: DW-370 Full-digital Trolley Ultrasound Scanner