




Infusion Pump

- Unique audible and visual alarm, safe and reliable, gain patents
- Compatible with any brands of iv sets after correct calibration
- Work-mode: ml/h and drop/min. Can switch freely



Features:

Audible and visual alarm for There are bubbles, Tube is blocked, Close the door, IV set Mismatched, Battery running out, Infusion finish, Equipment failure, Abnormally Operation, Software Exception, Low battery, forget the operations ect.
 Preset Solution Volume to greatly reduce the workload of nurses.
 Work mode: ml/h and drop/min can switch freely.
 Automatically record the settings of last infusion; BYS-820D can store 2000pcs of infusion records and it has a drug library.
 HD LCD Display, high capacity words, friendly user interface, dynamically display working status.
 OEM is available.

Specification	
Infusion Rate	ml/h mode: 1ml/h---1200ml/h(1ml/h step) drop/min mode: 1-266drop/min(1 drop/min step)
Flow Rate Accuracy	ml/h mode: Within±5%(after correct calibration) drop/min mode: Within±3%(after correct calibration)
Mechanical Precision	Within±2%
Purge Rate	100ml/h---1200ml/h(can be calibrated)
Volume Limit	1ml---9999ml
Total Infusion Volume	0---9999.9ml
Alarm Pressure	High: 800mmHg±200mmHg(106.7kPa±26.7kPa) Medium: 500mmHg±100mmHg(66.7kPa±13.3kPa) Low: 300mmHg±100mmHg(40.7kPa±13.3kPa)
Power Source	100V---240V 50/60Hz; Internal rechargeable Li battery, DC 12V
KVO Speed	1ml/h---5ml/h(1ml/h step)
Power Consumption	35VA
Infusion Set	Compatible with any brands of infusion sets after correct calibration
Dimension	140mm(L)×157mm(W)×220mm(H)
Net Weight	1.8Kg
Forbidden	Not support blood infusion
Package Reference	1pc Carton, 2.8Kg; 30cm(L)×21cm(W)×26cm(H) 2pcs Carton, 5.6Kg; 45cm(L)×25cm(W)×33cm(H) 3pcs Carton, 8.4Kg; 67cm(L)×25cm(W)×33cm(H)